



View

View allows providers to review imaging studies without ever leaving the EHR.

PRODUCT OVERVIEW

View is a SMART on FHIR app that opens right from the patient's chart with no separate portal, no second login, and no retyping information the EHR already has. From a single launch, a provider can review a patient's imaging studies in a non-diagnostic viewer, all without leaving the chart.

How View Works

View uses the SMART App Launch standard, the same framework EHR vendors use for their own embedded apps.

When a provider launches View:

1

The EHR passes launch context to Medicom, identifying the patient, the practitioner, and (optionally) any specific imaging studies.

2

Medicom completes the OAuth handshake with the EHR, confirming the user and organization.

3

View opens with a secure session to the Medicom platform, scoped to that patient and that provider.

Provider and Operational Benefits

Keeping providers in the EHR pays off well beyond convenience. It saves time, reduces risk, and speeds care.

FOR PROVIDERS

A provider opens View already authenticated and in the right patient's chart, so reviewing outside imaging takes one click instead of a separate login and a fresh search.

FOR OPERATIONAL TEAMS

There's no new credential system to manage and no desktop app to deploy. The EHR handles identity, and Medicom handles the imaging, so support tickets stay low.

In-Chart Clinical Power

Integrating external imaging management directly into the patient's chart allows clinical teams to process and review incoming studies without switching applications or interrupting care. Care teams can complete essential imaging tasks without leaving the clinical view.



Review & Reconcile With Checkpoint

Incoming outside studies land in a worklist next to their reports, so staff can link a study to an open order, create an order and store it, or store it without an order, based on how your organization sets policy.



Import From Physical Media With Data Import

Discs, USB drives, and folders import right from the patient's chart, so studies get checked against the right patient and can be matched to an open order from the start.



View Images Without Leaving the Chart

An embedded viewer opens studies alongside the worklist or in a separate tab, so no one exports a file or switches applications to see an image.

View runs inside the EHR your team already uses. [Schedule a customized walkthrough today.](#)

